



Chemical, Biological, Radiological & Nuclear Collective Protection

CBRN "Single Skin" Col Pro



Dual Use

UTS Systems "Single Skin" design eliminates the need for an additional CBRN Col Pro liner, increasing ease of set-up and decreasing logistics

- Manufactured with BondBarrier™ a Mil Spec 44103 fabric
- Utilizes Saint-Gobain Coretech®, a CBRN resistant barrier



Airlocks

Multi-Person Entrance (MPE)

- Double Bump Through Doors
- Accommodates two NATO litters with attendants or 12 people
- Purge time (99.9%) of less than 4 minutes with recirculation filter



Testing

- Chemical permeability live agent 72 hour testing against HD, GB, GD and Vx
- Chamber testing successfully completed including:
 - Over Pressure, Leakage, & Complexing
 - Airlock Purge (3-log reduction in concentration)
 - Static Vapor Challenge Test – Simulant: TMP (*Trimethyl Phosphate simulant for Sarin*)
 - Static Particulate Challenge Test (SPCT) – Simulant: BG (*Bacillus subtilis simulant for Anthrax*)
- Successful concept demonstration with Air Force Medical Service at Langley AFB



850.226.4301 • sales@uts-systems.com • www.uts-systems.com
36 Tupelo Ave SE • Ft. Walton Beach, Florida 32548