

CBRN Mass Personnel & Patient (MP&P) Decontamination (DECON)

Staging Area, DECON & Anteroom



- "Single Skin" Chem Bio resistant barrier
- Improved throughput using filtered force air flow and specialized flooring
- Compatible with Col Pro or standalone
- Shelter sleeve drop for quick removal of contaminated garments from interior

Integrated

- Containment berm integral to floor
- Electrical & Hoses housed in vertical/horizontal raceways
- Coiled sprayers for mixed solution & rinse and overhead spray bars
- Specially designed floor risers allow flow to sump pump & reduce slippage





Litter Conveyor

- UTS MP&P provides smooth processing for Non-Ambulatory Patients
- Secure locking, adjustable legs and leveling indicator
- Anodized aluminum rollers, nylon rod with stainless steel axle nuts-RUST FREE



DECON Water Heater

- UTS MP&P includes containerized dual on-demand diesel water heaters
- Allows for water to be drawn from:
 - Self contained 3,000 gallon tank
 - Pressurized feed from fire hydrant
 - Open source through inlet strainer
- Variable solution injection system

Waste Sump Pump

- Robust run dry sump powered by container
- Suction head reduces required space in shelter through one of four entry sleeves
- Ships securely inside container along with all hoses and associated equipment



